



## **ATESE ET-300 C OPERATION TABLE**

**With Zeroing Option**



+90 312 511 56 00



info@labcore.com.tr



Turgut Özal Mah. 2164 Sk.  
A Blok No:3 A-134  
Yenimahalle / Ankara

TABLE POSITIONS	
Head Section (Manual)	55'UP 90'down, Detachable
Lateral Tilt (Control)	25'Right 25' Left Lateral-Tilt position
Trendelenburg (Control)	30
Reverse Trendelenburg (Control)	25
Back Section (Control)	60°UP/40° Down (x15°)
Leg Section (Control)	30 Up . 90' Down.90' Side Sliding Detachable
Type	Electro-Mechanical(Control)
Slide (Optional) (Control)	300mm

TECHNICAL FEATURES	
MEMORY FOR POSITIONS	
ZERO POSITIONING	
BATTERY INDICATOR	
Minimum Height	710 mm x 1010 mm
Maximum Height	760 mm x 1060 mm
Tabletop Length	600 mm x 2150 mm
Head Section	225 mm x 355 mm
Back Section	530 mm x 625 mm
Sitting Base Section	530 mm x 530 mm
Leg Section (Right) (Piston)	255 mm x 710 mm
Leg Section (Left) (Piston)	255 mm x 710 mm
Power Input	230V 50Hz
Battery Charge Voltage	24V DC
Remote Control Voltage	5V DC
Engine Battery Voltage	220V AC, %a- 24V DC 10A
Battery	Dry Type / Maintenance Free
Battery Capacity	7Ah, 24V DC
Carrying Capacity	300 Kg
Battery Surviving Time	45 Min

+90 312 511 56 00

info@labcore.com.tr

Turgut Özal Mah. 2164 Sk.  
A Blok No:3 A-134  
Yenimahalle / Ankara