



TRACTION SET FOR OPERATING TABLE



+90 312 511 56 00



info@labcore.com.tr



Turgut Özal Mah. 2164 Sk.
A Blok No:3 A-134
Yenimahalle / Ankara

Model No	TS-I	Lifting Range	350mm
Traction Range	190mm	Material	304 Stainless Steel
Application Type	Any Operation Table	Warranty	1 year

Specifications

Height:800mm Lifting range:350mm Traction range:190mm

Main Function :

TSI is used for orthopedic traction, made of premium 304 stainless steel can use with our any operating table models (universal used)

Detailed Description:

- * Multifunctional universal joint
- *Space-saving storage of all parts
- * For all kind of fractures of the lower limbs
- *Easy transport and assembly thanks to trolley
- *Comfortable and precise positioning of the patient
- *Quick conversion of normal table to extension table

Technical Data:

Type	1006
Height(mm)	<=800
Lifting range (mm)	>=350
Tracation range	>=190
Tracation angle: rotation (up and down)	120
Horizontal rotation	360
Horizontal evagination of body frame	>=180
Extendable adjustment(back & forth)(mm)	>=400
Liftitng adjustment range of side (mm)	>=200

+90 312 511 56 00

info@labcore.com.tr

Turgut Özal Mah. 2164 Sk.
A Blok No:3 A-134
Yenimahalle / Ankara